

## US EPA Method 624 using the Teledyne Tekmar Atomx XYZ and Thermo Scientific™ ISQ™ 7000 Mass Spectrometry (MS) System Coupled with a Thermo Scientific™ TRACE™ 1310 Gas Chromatograph (GC)

Amy Nutter, Applications Chemist; Teledyne Tekmar

Page | 1

### Abstract

US EPA Method 624 was used to determine the concentration of volatile organic compounds (VOCs) in wastewater samples. The Teledyne Tekmar Atomx XYZ purge and trap (P&T) VOC sample preparation system combined with a Thermo Scientific ISQ 7000 Mass Spectrometry (MS) system coupled with a Thermo Scientific TRACE 1310 Gas Chromatograph (GC) was used to create a working calibration curve, method detection limits (MDLs) and initial demonstration of capability (IDC) accuracy and precision for target compounds. This study will demonstrate the ability of the Atomx XYZ's innovative moisture control system (MCS) to remove water vapor transferred to the GC/MS.

### Introduction

The Atomx XYZ is Teledyne Tekmar's most advanced P&T system and is based on the time-tested Atomx instrument platform. The concentrator's efficient trap cooling design reduces sample cycle time by as much as 14% over the previous model. Combined with its 84-position soil and water autosampler, the result is more samples tested per 12-hour period. The redesigned MCS improves water vapor removal by as much as 60%, thereby reducing peak interference and increasing GC column lifespan. In addition to other refinements, the Atomx XYZ incorporates a precision-machined valve manifold block to reduce potential leak sources and ensure the system is both reliable and robust.

### Sample Preparation

A working 50 ppm calibration standard was prepared in methanol from Restek® standards: 624 Calibration Mix #1 and Volatiles MegaMix™ and EPA Method 624. In total, the standard contained 31 compounds.

The water calibration curve was prepared from 0.5 ppb to 200 parts per billion (ppb) for all compounds. The relative response factor (RF) was calculated for each compound using one of the three internal standards: Bromochloromethane, 2-Bromo-1-chloropropane and 1,4-Dichlorobutane. Surrogate standards consisted of: Pentafluorobenzene, Fluorobenzene and 4-Bromofluorobenzene. Internal and surrogate standards were prepared together in methanol from Restek® standards at a concentration of 30 parts per million (ppm), after which 5 µL was then mixed with each 5 mL sample for a resulting concentration of 30 ppb.

Seven 0.5 ppb standards were prepared for MDLs and precision calculations. Also, seven 5 ppb standards were prepared for the IDC precision and accuracy calculations. All calibration, MDL and IDC samples were analyzed with the Atomx XYZ conditions in [Table I](#) and the GC/MS conditions in [Table II](#).

## Experimental Instrument Conditions

| Table I Teledyne Tekmar Atomx XYZ Water Method Conditions |                 |                                 |                 |
|---|-----------------|---------------------------------|-----------------|
| Purge   | Variable        | Desorb                          | Variable        |
| Valve Oven Temp   | 140 °C          | Methanol Needle Rinse           | Off             |
| Transfer Line Temp  | 140 °C          | Methanol Needle Rinse Volume    | 0.00 mL         |
| Sample Mount Temp   | 90 °C           | Water Needle Rinse Volume       | 7.00 mL         |
| Water Heater Temp   | 90 °C           | Sweep Needle Time               | 0.25 min        |
| Sample Vial Temp  | 20 °C           | Dry Purge Temp                  | 20 °C           |
| Soil Valve Temp   | 100 °C          | Desorb Preheat Temp             | 245 °C          |
| Standby Flow  | 10 mL/min       | GC Start Signal                 | Begin Desorb    |
| Condensate Ready Temp                                     | 45 °C           | Desorb Time                     | 2.00 min        |
| Purge Ready Temp  | 40 °C           | Drain Flow                      | 300 mL/min      |
| <b>Purge</b>  | <b>Variable</b> | Desorb Temp                     | 250 °C          |
| Sample Equilibrate Time                                   | 0.00 min        | <b>Bake</b>                     | <b>Variable</b> |
| Pre-sweep Time  | 0.25 min        | Methanol Glass Rinse            | Off             |
| Prime Sample Fill Volume                                  | 3.00 mL         | Number of Methanol Glass Rinses | 0               |
| Sample Volume   | 5.00 mL         | Methanol Glass Rinse Volume     | 0.00 mL         |
| Sweep Sample Time   | 0.25 min        | Water Bake Rinses               | 1               |
| Sweep Sample Flow   | 100 mL/min      | Water Bake Rinse Volume         | 7.00 mL         |
| Spurge Vessel Heater                                      | Off             | Bake Rinse Sweep Time           | 0.25 min        |
| Spurge Vessel Temp  | 20 °C           | Bake Rinse Sweep Flow           | 100 mL/min      |
| Pre-purge Time  | 0.00 min        | Bake Rinse Drain Time           | 0.40 min        |
| Pre-purge Flow  | 0 mL/min        | Bake Time                       | 2.00 min        |
| Purge Time  | 11.00 min       | Bake Flow                       | 200 mL/min      |
| Purge Flow  | 40 mL/min       | Bake Temp                       | 260 °C          |
| Purge Temp  | 20 °C           | Condensate Bake Temp            | 200 °C          |
| Condensate Purge Temp                                     | 20 °C           |                                 |                 |
| Dry Purge Time  | 0.00 min        | <b>Trap</b>                     | #9              |
| Dry Purge Flow  | 100 mL/min      | <b>Purge Gas</b>                | Nitrogen        |

| Table II Thermo Scientific TRACE 1310 GC and ISQ 7000 MS System Conditions |   |
|--|---|
| Thermo Scientific TRACE 1310 GC Conditions                                 |   |
| Column   | TG VMS, 20 m x 0.18 mm, 1µm Film, Helium – 0.8 mL/min                                     |
| Oven Profile   | 35 °C, 3 min, 12 °C/min to 85 °C, 25 °C/min to 225 °C, 2 min Hold, Run Time 14.767 min    |
| Inlet  | 200 °C, 50:1 Split, Purge Flow 0.5 mL/min   |
| Thermo Scientific ISQ 7000 MS Conditions                                   |   |
| Temp   | Transfer Line 230 °C; Ion Source 280 °C   |
| Scan   | Range 35 <i>amu</i> to 260 <i>amu</i> , Solvent Delay 0.50 min, Dwell/Scan Time 0.15 sec. |
| Current  | Emission Current 25 µA, Gain 3.00E+005  |

## Results


The relative standard deviation (%RSD) of the RFs for the calibration curve, MDL, precision and IDC data are shown in Table III. Figure 1 displays a 5 ppb standard, indicating excellent peak resolution with minimal water inference for VOCs standards.

| Table III US EPA Method 8260 Water Calibration, Accuracy and Precision Data |                |                         |            |                                       |                             |                  |  |                      |
|---|----------------|-------------------------|------------|---------------------------------------|-----------------------------|------------------|--|----------------------|
| Compound  | Calibration    |                         |            | Method Detection Limit (n=7, 0.5 ppb) |                             |                  | Initial Demonstration of Capability (n=7, 5 ppb) |                      |
|   | Retention Time | Linearity RF (≤20% RSD) | Average RF | MDL (ppb)                             | Average Concentration (ppb) | Precision (≤20%) | Accuracy   | Precision (≤20% RSD) |
| Chloromethane   | 1.51           | 17.8                    | 1.36       | 0.31                                  | 0.54                        | 18.4             | 83   | 13.2                 |
| Vinyl Chloride  | 1.58           | 5.61                    | 0.672      | 0.12                                  | 0.43                        | 8.91             | 91   | 12.9                 |
| Bromomethane  | 1.85           | 16.7                    | 0.349      | 0.17                                  | 0.69                        | 7.61             | 101  | 11.3                 |
| Chloroethane <sup>1</sup>   | 1.96           | 16.2                    | 0.382      | 0.28                                  | 0.68                        | 12.9             | 119  | 12.5                 |
| Trichlorofluoromethane  | 2.09           | 5.76                    | 0.610      | 0.15                                  | 0.44                        | 10.9             | 90   | 13.4                 |
| 1,1-Dichloroethene  | 2.52           | 6.09                    | 0.396      | 0.11                                  | 0.44                        | 8.01             | 94   | 14.0                 |
| Methylene Chloride  | 3.07           | 14.5                    | 0.576      | 0.09                                  | 0.56                        | 5.03             | 95   | 13.1                 |
| trans-1,2-Dichloroethene  | 3.22           | 5.83                    | 0.449      | 0.12                                  | 0.50                        | 7.26             | 97   | 13.8                 |
| 1,1-Dichloroethane  | 3.84           | 5.41                    | 1.09       | 0.09                                  | 0.45                        | 6.13             | 95   | 12.9                 |
| Bromochloromethane (IS)   | 4.57           |                         |            |                                       |                             |                  |  |                      |
| Chloroform  | 4.67           | 6.42                    | 0.793      | 0.08                                  | 0.44                        | 5.97             | 96   | 12.8                 |
| Carbon Tetrachloride  | 4.76           | 13.7                    | 1.66       | 0.09                                  | 0.38                        | 7.71             | 101  | 14.3                 |
| 1,1,1-Trichloroethane   | 4.83           | 8.67                    | 2.17       | 0.11                                  | 0.41                        | 8.25             | 100  | 13.8                 |
| Benzene   | 5.20           | 6.55                    | 7.24       | 0.09                                  | 0.44                        | 6.73             | 94   | 12.8                 |
| Pentafluorobenzene (SURR)   | 5.34           | 8.64                    | 3.78       |                                       | 30.0                        | 2.19             | 99   | 2.49                 |
| 1,2-Dichloroethane  | 5.41           | 13.9                    | 0.241      | 0.29                                  | 0.60                        | 15.4             | 105  | 11.9                 |
| Fluorobenzene (SURR)  | 5.62           | 1.75                    | 7.39       |                                       | 29.8                        | 0.994            | 100  | 0.726                |
| Trichloroethene   | 5.78           | 10.2                    | 1.36       | 0.12                                  | 0.44                        | 8.49             | 100  | 14.6                 |
| 1,2-Dichloropropane <sup>3</sup>  | 6.29           | 10.1                    | 0.233      | 0.32                                  | 0.94                        | 10.7             | 94   | 16.1                 |
| Bromodichloromethane  | 6.38           | 9.20                    | 0.154      | 0.18                                  | 0.42                        | 13.7             | 100  | 12.8                 |
| 2-Chloroethyl Vinyl Ether   | 6.99           | 3.90                    | 0.234      | 0.17                                  | 0.46                        | 11.8             | 91   | 13.9                 |
| trans-1,3-Dichloropropene   | 7.02           | 10.4                    | 2.73       | 0.10                                  | 0.40                        | 7.75             | 88   | 14.0                 |

**Table III US EPA Method 8260 Water Calibration, Accuracy and Precision Data**

| Compound                               | Calibration       |                            |               | Method Detection Limit<br>(n=7, 0.5 ppb) |                                |                     | Initial Demonstration of<br>Capability (n=7, 5 ppb) |                         |
|--|-------------------|----------------------------|---------------|--|--------------------------------|---------------------|---|-------------------------|
|  | Retention<br>Time | Linearity RF<br>(≤20% RSD) | Average<br>RF | MDL<br>(ppb)                             | Average<br>Concentration (ppb) | Precision<br>(≤20%) | Accuracy  | Precision<br>(≤20% RSD) |
| Toluene                                | 7.24              | 7.18                       | 3.81          | 0.12                                     | 0.45                           | 8.42                | 94  | 13.4                    |
| Tetrachloroethene                      | 7.59              | 6.07                       | 0.315         | 0.13                                     | 0.44                           | 9.67                | 87  | 13.3                    |
| 2-Bromo-1-chloropropane (IS)           | 7.64              |                            |               |  |                                |                     |   |                         |
| cis-1,3-Dichloropropene <sup>2,3</sup> | 7.65              | 17.7                       | 0.721         | 0.22                                     | 1.2                            | 5.72                | 81  | 11.8                    |
| 1,1,2-Trichloroethane                  | 7.80              | 9.54                       | 0.345         | 0.07                                     | 0.43                           | 4.94                | 91  | 14.0                    |
| Dibromochloromethane                   | 7.93              | 12.2                       | 0.276         | 0.09                                     | 0.41                           | 6.71                | 87  | 13.6                    |
| Chlorobenzene                          | 8.55              | 8.35                       | 1.09          | 0.10                                     | 0.45                           | 6.99                | 88  | 12.3                    |
| Ethylbenzene                           | 8.58              | 6.55                       | 0.613         | 0.12                                     | 0.42                           | 9.14                | 84  | 13.9                    |
| Bromoform                              | 9.07              | 7.80                       | 0.259         | 0.07                                     | 0.41                           | 5.28                | 83  | 14.4                    |
| 4-Bromofluorobenzene (SURR)            | 9.45              | 4.37                       | 0.772         |  | 28.5                           | 2.19                | 95  | 2.18                    |
| 1,4-Dichlorobutane (IS)                | 9.57              |                            |               |  |                                |                     |   |                         |
| 1,1,2,2-Tetrachloroethane              | 9.61              | 9.60                       | 0.038         | 0.18                                     | 0.47                           | 11.9                | 88  | 14.3                    |
| 1,3-Dichlorobenzene                    | 10.17             | 10.9                       | 1.20          | 0.14                                     | 0.50                           | 9.01                | 87  | 13.3                    |
| 1,4-Dichlorobenzene                    | 10.24             | 10.7                       | 1.24          | 0.13                                     | 0.51                           | 8.10                | 86  | 12.8                    |
| 1,2-Dichlorobenzene                    | 10.52             | 10.0                       | 1.17          | 0.13                                     | 0.48                           | 8.47                | 89  | 13.3                    |

1. Calibration curve 2-200 ppb.
2. Calibration curve 1-200 ppb.
3. MDL with 1 ppb.

**Figure 1** Total Ion Chromatogram of a Water Method 5 ppb VOC Standard Indicating Consistent Peak Shapes for all Compounds with Minimal Water Interference.


## Conclusion

This study demonstrates the capability of the Teledyne Tekmar Atomx XYZ P&T system to process VOCs in wastewater samples following US EPA Method 624 with detection by a Thermo Scientific ISQ 7000 Mass Spectrometry (MS) system coupled with a Thermo Scientific TRACE 1310 Gas Chromatograph (GC). The %RSD of the calibration curve passed all method requirements with no interference from excessive water. Furthermore, MDL, precision and IDC for seven 0.5 ppb and seven 5 ppb standards showed no interference from excessive water.

By making additional, appropriate changes to the GC oven temperature program, the GC/MS cycle time may also be reduced, increasing laboratory throughput in a 12-hour period.

## References

1. Appendix A to Part 136 - Methods for Organic Chemical Analysis of Municipal and Industrial Wastewater - Method 624: Purgeables; US EPA, Promulgated 1984. [Online] [https://www.epa.gov/sites/production/files/2015-10/documents/method\\_624\\_1984.pdf](https://www.epa.gov/sites/production/files/2015-10/documents/method_624_1984.pdf) (accessed March 26, 2020).