

---

# Conducting Batch Reactions at Relevant Geologic Carbon Sequestration Conditions

---



## Syringe Pump Application Note AN42

### Overview


Geologic carbon sequestration (GCS) is needed to mitigate global warming caused by increased atmospheric CO<sub>2</sub> concentrations. For safe and efficient GCS operations the geochemical reactions with CO<sub>2</sub>-brine-rock under relevant GCS conditions needs to be studied. For this application CO<sub>2</sub>-brine-biotite interactions were studied and fibrous illite formation was observed after reaction for only 3 hours. The fibrous illite detached and released into solution after longer reaction time, indicating potential permeability decrease in the aquifer, as fibrous illite can clog the pore spaces affecting further CO<sub>2</sub> injections.

### Experiment

At GCS sites, there is high temperature and high pressure, with high salt concentrations (especially NaCl) in the brine. To mimic the GCS conditions, as shown in the Figure 1, gas (e.g, CO<sub>2</sub>) can be injected into the high temperature reactor by the Teledyne ISCO syringe pump. This single pump set-up was used for its ability to maintain a constant pressure. During reaction, CO<sub>2</sub> may dissolve into the solution. However, it is important to leave the valve between the syringe pump and the reactor open and set the syringe pump to be in constant pressure mode.

During reaction, CO<sub>2</sub> in the solution can be consumed and thus CO<sub>2</sub> can be continuously pumped into the reactor to maintain the set pressure. After reaction for different elapsed time (e.g., 3, 5, 8, 17, 22, 44, 70 hrs and 6 days), the valve connecting the syringe pump and the reactor was closed to prevent further CO<sub>2</sub> injection into the reactor. Then the gas release valve of the reactor was open to release the CO<sub>2</sub>. Meanwhile, the reactor was cooled down to room temperature.

After releasing the pressure and lower the temperature, the reactor was opened, the dissolved aqueous ion concentrations was analyzed by ICP-MS, to represent the biotite dissolution. The surfaces of the reacted biotite was analyzed by AFM, SEM, to study the surface morphology changes and secondary mineral formations.



rock/water ratio: 1 piece of biotite flake and 4 ml solution, 0.015 g/4 ml

Elapsed time: 3,5,8,17,22,44,70,96,120,144 h.


**Figure 1: Single pump used to mimic GCS conditions**

#### **Note**


Experiments conducted in DI water were serving as a control. For clarification, all images shown in Figure 2 were collected after reaction in 1 M NaCl solutions at 95° C and 102 atm CO<sub>2</sub>. Experiments were also conducted at lower temperatures. At 35° C the kinetics of illite formations slowed down, but illites formed at all temperatures.

## Conclusions

As shown in Figures 2 and 3, the formation and mobilization of fibrous illite from biotite under GCS conditions can lower the permeability of the aquifer. So it may cause problems during CO<sub>2</sub> injection.



**Figure 2: AFM observation of fibrous illite formation on reacted basal surface**



**Detachment of fibrous illite and biotite flakes from cracked surfaces**

February 13, 2014; revised November 8, 2023

#### **Teledyne ISCO**

P.O. Box 82531, Lincoln, Nebraska, 68501 USA  
Toll-free: (800) 228-4373 • Phone: (402) 464-0231 • Fax: (402) 465-3091  
www.teledyneisco.com

Teledyne ISCO is continually improving its products and reserves the right to change product specifications, replacement parts, schematics, and instructions without notice.

EAR99 Technology Subject to Restrictions Contained on the Cover Page

