

Advanced Topics for Use of RediSep[®] Specialty Media Columns in Flash Chromatography



Veronica Thomason

Chromatography Applications Chemist
E-mail: vthomason@teledyne.com

Teledyne Isco, Inc.

P.O. Box 82531
Lincoln, NE 68501

Web: www.isco.com

Abstract

RediSep[®] specialty media columns offer synthetic organic chemists practical tools for obtaining optimal purification by flash chromatography. RediSep Cyano columns are featured as an efficient alternative to normal phase silica gel for high pK_a compound separations. In addition, a practical approach for using specialty media such as RediSep Cyano and RediSep Amine columns under reversed phase conditions is discussed.






Figure 1.
Solvents demonstrate different strengths with varying media polarity

Separation of High pK_a Compounds on RediSep Cyano

Silica gel is a widely used stationary phase for normal phase purification of compounds with low to medium polarity. Yet, inefficiencies of silica gel are often encountered by the synthetic chemist when purifying basic or medium-to-high polarity compounds. High pK_a compounds such as 3-(dimethylaminomethyl) indole (**1**) and 2-(dimethylaminomethyl)-3-hydroxypyridine (**2**), are often fully retained on silica gel due to the presence of acidic silanol groups, requiring the use of basic modifiers, such as ammonia, ammonium hydroxide, or triethylamine in the mobile phase. In addition, such relatively polar compounds have a strong affinity to the polar silanol groups and require elution through silica gel with a mobile phase of higher solvent strength than ethyl acetate in hexane, such as 10% methanol in dichloromethane (Figure 3, top).



(1)



(2)

RediSep Cyano functionalized silica is a neutral media with a short alkyl chain end capped with a weakly polar nitrile group. High pK_a compounds such as 3-(dimethylaminomethyl) indole (**1**) and 2-(dimethylaminomethyl)-3-hydroxypyridine (**2**) are not retained as strongly (Figure 3, bottom), thereby eliminating the need for basic modifiers and the use of chlorinated solvent systems.

Alternatively, RediSep Cyano columns also can be used as an efficient alternative to C-18 silica for the separation of compounds of medium to high polarity. The short alkyl chain does not retain moderately polar compounds as strongly in reversed phase mode like C-18 silica with its long alkyl chain. Strong retention of less polar compounds on C-18 can result in broader peak shape and longer elution time. These effects are seen with the broader peak shape of the hydroxypyridine derivative on C-18 and the indole derivative not eluting within the designated time frame (Figure 4, bottom). Separation of moderately polar compounds is enhanced on RediSep Cyano due to the decrease of retentive properties while maintaining good selectivity.

Table 1 Experimental



Instrument:	Teledyne Isco's CombiFlash® Rf
Columns:	RediSep Flash 4g (Silica gel), RediSep Cyano 5g, RediSep C-18 Reversed Phase 4.3g
Sample mass:	40mg
Flow rate:	16ml/min
Detection wavelength:	280nm



Figure 2.
Teledyne Isco's CombiFlash Rf
automated Flash chromatography system



RediSep Silica Gel (highly polar)

Other key property:
Slightly acidic



RediSep Cyano (weakly polar)

Other key properties:
▪ Non-acidic
▪ Short alkyl chain



A: CH₂Cl₂
B: CH₃OH

Normal
Phase



A: Hexane
B: EtOAc

Figure 3.
Normal phase separation of high pK_a compounds

RediSep Cyano (weakly polar)

Other key properties:

- Non-acidic
- Short alkyl chain





A: H₂O
B: CH₃CN

Reversed
Phase

RediSep C-18 (non-polar)

Other key property:
Long alkyl chain



A: H₂O
B: CH₃CN


Figure 4.
Reversed phase separation of high pK_a compounds

Separation of Xanthine Derivatives on RediSep Amine


Synthetic organic chemists must develop unique separation methods for compounds that are difficult to separate. Developing a practical separation method requires the appropriate stationary phase and mobile phase combination to create the selectivity needed for an effective separation of such a sample. It is general practice that in reversed phase separation, the non-polar stationary phase is used with the mobile phase increasing in strength with decreasing solvent polarity (*e.g.*, water to acetonitrile), thus eluting the most polar compounds first. As illustrated with the separation of xanthine derivatives on C-18 silica (Figure 5, top), the solvent strength increased from water to acetonitrile and resulted in a partial separation with the relatively more polar compound, theophylline (**4**), eluting first. Conversely, for the separation of the xanthine derivatives using the RediSep Amine stationary phase, it can be seen that the compounds have low retention when the same gradient profile is used, eluting the relatively less polar compound, caffeine (**3**), first (Figure 5, bottom). Such a drastic change in result would suggest a review of the selectivity of the mobile phase with the amino stationary phase for the separation of xanthine derivatives.

The RediSep Amine stationary phase is functionalized silica with a short alkyl chain end capped with a polar amino group. In an attempt to use the RediSep Amine column in a reverse-phase separation mode (*i.e.*, water as the weak solvent), a fast elution of the xanthine derivatives resulted with an inversion of compound elution order. It is general practice that in normal phase separation, solvent strength increases with increasing solvent polarity, typically eluting the less polar compounds first. Using these reversed phase solvents on the

RediSep Amine stationary phase resulted in complete separation of the xanthine derivatives when using the appropriate solvent, acetonitrile, as the weak solvent; eluting the relatively less polar compound, caffeine (**3**), first as would be expected for a normal phase mode of separation. It can be seen that retention of the xanthine derivatives was increased when beginning the gradient profile with the appropriate “weak” solvent, acetonitrile.



(3)





(4)

Table 2 Experimental

Instrument:	Teledyne Isco's CombiFlash® Rf
Columns:	RediSep C-18 Reversed Phase 4.3g, RediSep Amine 4.7g
Sample mass:	40mg
Flow rate:	16ml/min
Detection wavelength:	280nm

RediSep C-18
(non-polar)

Other key property:
Long alkyl chain





H₂O = weak

A: H₂O
B: CH₃CN

RediSep Amine
(slightly polar)

Other key properties:
▪ Basic
▪ Short alkyl chain





H₂O = strong

Figure 5.
Comparison of separations on RediSep C-18 and
Amine columns with H₂O and CH₃CN mobile phase

RediSep Amine
(slightly polar)

Other key properties:

- Basic
- Short alkyl chain



CH₃CN = weak

A: H₂O
B: CH₃CN

Figure 6.
Normal phase separation with H₂O and CH₃CN
on RediSep Amine column. Gradient profile
initiated with the weak solvent CH₃CN.



Figure 7.
Teledyne Isco's RediSep C-18 reversed phase
chromatography columns



Summary

It is demonstrated that the RediSep Cyano stationary phase is an effective alternative to silica for the purification of basic compounds due to the media being far less acidic than normal phase silica. In addition, RediSep Cyano can be used as an efficient alternative to C-18 silica for the separation of compounds of medium polarity by decreasing elution time and enhancing peak shape.

It has been illustrated that RediSep Amine stationary phase is favorable for performing separations in a normal phase mode employing a wide variety of solvents—including those typically used in reversed phase separations.

© 2007 Teledyne Isco, Inc.

RediSep® and CombiFlash® are registered trademarks of Teledyne Isco, Inc. All other brand and product names are trademarks or registered trademarks of their respective holders.



TELEDYNE ISCO

A Teledyne Technologies Company

P.O. Box 82531, Lincoln, Nebraska, 68501 USA

Toll-free: (800) 228-4373 ■ Phone: (402) 464-0231 ■ Fax: (402) 465-3091

E-mail: IscoInfo@teledyne.com